




# FWxC<sup>®</sup>

## Motorized Filter Wheels

### SDK Manual





# Table of Contents


<b>Chapter 1</b>	<b>Introduction .....</b>	<b>1</b>
<b>Chapter 2</b>	<b>C++ Software Development Kit .....</b>	<b>3</b>
<b>2.1.</b>	<b><i>FWxCCommand.h File Reference.....</i></b>	<b>3</b>
2.1.1.	Functions.....	3
2.1.2.	Function Documentation .....	3
<b>Chapter 3</b>	<b>Python Software Development Kit.....</b>	<b>11</b>
<b>3.1.</b>	<b><i>FWxC_COMMAND_LIB Namespace Reference .....</i></b>	<b>11</b>
3.1.1.	Functions.....	11
3.1.2.	Function Documentation .....	11
<b>Chapter 4</b>	<b>LabVIEW Software Development Kit .....</b>	<b>15</b>
<b>4.1.</b>	<b><i>How to install .....</i></b>	<b>15</b>
<b>4.2.</b>	<b><i>How to find VI.....</i></b>	<b>15</b>
<b>4.3.</b>	<b><i>How to use .....</i></b>	<b>16</b>

---


## Chapter 1 Introduction

User can start software development in C/C++ develop environment, LabVIEW and Python.


The software development interface can be found in the start menu.



or by clicking *Help* in software menu.



In this directory, you will find the support files for software development, as shown below.



## Chapter 2 C++ Software Development Kit

User can start software development with **FilterWheel102\_win32.dll** or **FilterWheel102\_win64.dll** in C/C++ development environment which can be found in **Thorlabs\_FWxC\_C++SDK** under \Document directory. The corresponding header file is also in **Thorlabs\_FWxC\_C++SDK** under \Document directory.

Copy **FilterWheel102\_win32.dll** or **FilterWheel102\_win64.dll** to your program folder, and make sure the library file and exe file are in the same folder.

Below is the description of the header file FWxCCommand.h

### 2.1. FWxCCommand.h File Reference

#### 2.1.1.Functions

- **DllExport int List** (unsigned char \*serialNo, unsigned int length)  
*list all the possible port on this computer.*
- **DllExport int Open** (char \*serialNo, int nBaud, int timeout)  
*open port function.*
- **DllExport int IsOpen** (char \*serialNo)  
*check opened status of port*
- **DllExport int Close** (int hdl)  
*close current opened port*
- **DllExport int Read** (int hdl, unsigned char \*b, int limit)
- **DllExport int Write** (int hdl, char \*b, int size)
- **DllExport int Set** (int hdl, char \*c, int size)
- **DllExport int Get** (int hdl, unsigned char \*c, unsigned char \*d)
- **DllExport int Purge** (int hdl, int flag)  
*Purge the RX and TX buffer on port.*
- **DllExport int SetPosition** (int hdl, int pos)
- **DllExport int SetPositionCount** (int hdl, int count)
- **DllExport int SetSpeedMode** (int hdl, int mode)
- **DllExport int SetTriggerMode** (int hdl, int mode)
- **DllExport int SetSensorMode** (int hdl, int mode)
- **DllExport int Save** (int hdl)
- **DllExport int GetPosition** (int hdl, int \*pos)
- **DllExport int GetPositionCount** (int hdl, int \*poscount)
- **DllExport int GetSpeedMode** (int hdl, int \*speed)
- **DllExport int GetTriggerMode** (int hdl, int \*triggermode)
- **DllExport int GetSensorMode** (int hdl, int \*sensormode)
- **DllExport int GetTimeToCurrentPos** (int hdl, int \*time)
- **DllExport int GetId** (int hdl, char \*d)

#### 2.1.2.Function Documentation

##### **DllExport int Close (int hdl)**

close current opened port

##### **Parameters**

<i>hdl</i>	handle of port.
------------	-----------------

**Returns**

0: success; negtive number : failed.

**DllExport int Get (int hdl, unsigned char \* c, unsigned char \* d)**

set command to device according to protocol in manual and get the return string.

make sure the port was opened SUCCESSful before call this function.

make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>c</i>	input command string (<255)
<i>d</i>	output string (<255)

**Returns**

0: SUCCESS; negative number: failed.

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

**DllExport int GetId (int hdl, char \* d)**

get the fiterwheel id

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>d</i>	output string (<255)

**Returns**

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

**DllExport int GetPosition (int hdl, int \* pos)**

get the fiterwheel current position

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>pos</i>	fiterwheel actual position

**Returns**

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

### ***DllExport int GetPositionCount (int hdl, int \* poscount)***

get the fiterwheel current position count

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

#### ***Parameters***

<i>hdl</i>	handle of port.
<i>poscount</i>	fiterwheel actual position count

#### ***Returns***

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

### ***DllExport int GetSensorMode (int hdl, int \* sensormode)***

get the fiterwheel current sensor mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

#### ***Parameters***

<i>hdl</i>	handle of port.
<i>sensormode</i>	fiterwheel actual sensor mode:0, Sensors turn off;1, Sensors remain active

#### ***Returns***

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

### ***DllExport int GetSpeedMode (int hdl, int \* speed)***

get the fiterwheel current speed mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

#### ***Parameters***

<i>hdl</i>	handle of port.
<i>speed</i>	fiterwheel actual speed mode:0,slow speed:1,high speed

#### ***Returns***

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;



**DllExport int GetTimeToCurrentPos (int hdl, int \* time)**

get the fiterwheel current sensor mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>time</i>	the time from last position to current position

**Returns**

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

**</param> DllExport int GetTriggerMode (int hdl, int \* triggermode)**

get the fiterwheel current trigger mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>triggermode</i>	fiterwheel actual trigger mode:tr0, input mode;1, output mode

**Returns**

0: success;

0xEA: CMD\_NOT\_DEFINED;

0xEB: time out;

0xED: invalid string buffer;

**DllExport int IsOpen (char \* serialNo)**

check opened status of port

**Parameters**

<i>serialNo</i>	serial number of the device to be checked.
-----------------	--

**Returns**

0: port is not opened; 1 : port is opened.

**DllExport int List (unsigned char \* serialNo, unsigned int length)**

list all the possible port on this computer.

**Parameters**

<i>serialNo</i>	port list returned string include serial number and device descriptor, separated by comma
<i>length</i>	max size of buf

**Returns**

non-negative number: number of device in the list; negative number: failed.

**DllExport int Open (char \* serialNo, int nBaud, int timeout)**

open port function.

**Parameters**

<i>serialNo</i>	serial number of the device to be opened, use GetPorts function to get exist list first.
<i>nBaud</i>	bit per second of port
<i>timeout</i>	set timeout value in (s)

**Returns**

non-negative number: hdl number returned successfully; negtive number : failed.

**DllExport int Purge (int hdl, int flag)**

Purge the RX and TX buffer on port.

**Parameters**

<i>hdl</i>	handle of port.
<i>flag</i>	

FT\_PURGE\_RX: 0x01

FT\_PURGE\_TX: 0x02

**Returns**

0: SUCCESS; negative number: failed.

**DllExport int Read (int hdl, unsigned char \* b, int limit)**

read string from device through opened port.

make sure the port was opened SUCCESSful before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>b</i>	returned string buffer
<i>limit</i>	

ABS(limit): max length value of b buffer.

SIGN(limit) == 1 : wait RX event until time out value expired;

SIGN(limit) == -1: INFINITE wait event untill RX has data;

**Returns**

non-negative number: size of actual read data in byte; negative number: failed.

**DllExport int Save (int hdl)**

save all the settings as default on power up

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
------------	-----------------

**Returns**

0: success;  
0xEA: CMD\_NOT\_DEFINED;  
0xEB: time out;  
0xED: invalid string buffer;

**DllExport int Set (int hdl, char \* c, int size)**

set command to device according to protocol in manual.  
make sure the port was opened SUCCESSful before call this function.  
make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>c</i>	input command string
<i>size</i>	lenth of input command string (<255)

**Returns**

0: SUCCESS; negative number: failed.  
0xEA: CMD\_NOT\_DEFINED;  
0xEB: time out;  
0xED: invalid string buffer;

**DllExport int SetPosition (int hdl, int pos)**

set fiterwheel's position  
make sure the port was opened successful before call this function.  
make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>pos</i>	fiterwheel position

**Returns**

0: success;  
0xEA: CMD\_NOT\_DEFINED;  
0xEB: time out;  
0xED: invalid string buffer;

**DllExport int SetPositionCount (int hdl, int count)**

set fiterwheel's position count  
make sure the port was opened successful before call this function.  
make sure this is the correct device by checking the ID string before call this function.

**Parameters**

<i>hdl</i>	handle of port.
<i>count</i>	fiterwheel PositionCount

**Returns**

0: success;

0xEA: CMD\_NOT\_DEFINED;  
 0xEB: time out;  
 0xED: invalid string buffer;

### ***DllExport int SetSensorMode (int hdl, int mode)***

set fiterwheel's sensor mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

#### ***Parameters***

<i>hdl</i>	handle of port.
<i>mode</i>	fiterwheel sensor mode:sensors=0 Sensors turn off when wheel is idle to eliminate stray light;sensors=1 Sensors remain active

#### ***Returns***

0: success;  
 0xEA: CMD\_NOT\_DEFINED;  
 0xEB: time out;  
 0xED: invalid string buffer;

### ***DllExport int SetSpeedMode (int hdl, int mode)***

set fiterwheel's speed mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

#### ***Parameters***

<i>hdl</i>	handle of port.
<i>mode</i>	fiterwheel speed mode:speed=0 Sets the move profile to slow speed:speed=1 Sets the move profile to high speed

#### ***Returns***

0: success;  
 0xEA: CMD\_NOT\_DEFINED;  
 0xEB: time out;  
 0xED: invalid string buffer;

### ***DllExport int SetTriggerMode (int hdl, int mode)***

set fiterwheel's trigger mode

make sure the port was opened successful before call this function.

make sure this is the correct device by checking the ID string before call this function.

#### ***Parameters***

<i>hdl</i>	handle of port.
<i>mode</i>	fiterwheel trigger mode:trig=0 Sets the external trigger to the input mode, Respond to an active low pulse by advancing position by 1;trig=1 Sets the external trigger to the output mode, Generate an active high pulse when selected position arrived at

#### ***Returns***

0: success;  
0xEA: CMD\_NOT\_DEFINED;  
0xEB: time out;  
0xED: invalid string buffer;

**DllExport int Write (int hdl, char \* b, int size)**

write string to device through opened port.  
make sure the port was opened SUCCESSful before call this function.

**Parameters**






<i>hdl</i>	handle of port.
<i>b</i>	input string
<i>size</i>	size of string to be written.

**Returns**

non-negative number: number of bytes written; negative number: failed..

## Chapter 3 Python Software Development Kit

Python 3.6 or above is required. User can import **FWxC\_COMMAND\_LIB.py** to your python project, that's the wrapper for **FWxC\_COMMAND\_LIB(FilterWheel102\_win32.dll** in C/C++ development environment ). Copy **FilterWheel102\_win32.dll** to your program folder, and make sure the library file and **FWxC\_COMMAND\_LIB.py** file are in the same folder. The “ **FWxC\_COMMAND\_LIB\_EXAMPLE.py**” is example code for how to use the python APIs.

 <b>_pycache_</b>	12/4/2020 4:20 PM	File folder	
 <b>FilterWheel102_win32.dll</b>	11/11/2020 2:48 PM	Application extens	223 KB
 <b>FWxC_COMMAND_LIB.py</b>	12/4/2020 4:20 PM	Python source file	9 KB
 <b>FWxC_COMMAND_LIB_EXAMPLE.py</b>	12/4/2020 4:33 PM	Python source file	4 KB
 <b>Thorlabs_FWxC_PythonSDK.pyproj</b>	12/4/2020 4:14 PM	Python Project	2 KB

User can also replace the reference win32 lib to x64 lib for 64-bit application.

### 3.1. FWxC\_COMMAND\_LIB Namespace Reference

#### 3.1.1.Functions

- def **FWxCListDevices** ()
- def **FWxCOpen** (serialNo, nBaud, timeout)
- def **FWxCIsOpen** (serialNo)
- def **FWxCClose** (hdl)
- def **FWxCSetPosition** (hdl, pos)
- def **FWxCSetPositionCount** (hdl, count)
- def **FWxCSetSpeedMode** (hdl, spmode)
- def **FWxCSetTriggerMode** (hdl, trimode)
- def **FWxCSetSensorMode** (hdl, senmode)
- def **FWxCSave** (hdl)
- def **FWxCGetId** (hdl, value)
- def **FWxCGetPosition** (hdl, pos)
- def **FWxCGetPositionCount** (hdl, poscount)
- def **FWxCGetSpeedMode** (hdl, spemode)
- def **FWxCGetTriggerMode** (hdl, triggermode)
- def **FWxCGetSensorMode** (hdl, sensormode)

#### 3.1.2.Function Documentation

**def FWxC\_COMMAND\_LIB.FWxCClose ( hdl)**

```
Close opened FWxC device
Args:
    hdl: the handle of opened FWxC device
Returns:
    0: Success; negative number: failed.
```

**def FWxC\_COMMAND\_LIB.FWxCGetId ( hdl, value)**

```
get the FWxC id
Args:
    hdl: the handle of opened FWxC device
    value: the model number, hardware and firmware versions
Returns:
```

---

```
0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

### **def FWxC\_COMMAND\_LIB.FWxCGetPosition ( hdl, pos)**

```
get the fiterwheel current position
Args:
    hdl: the handle of opened FWxC device
    pos: fiterwheel actual position
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

### **def FWxC\_COMMAND\_LIB.FWxCGetPositionCount ( hdl, poscount)**

```
get the fiterwheel current position count
Args:
    hdl: the handle of opened FWxC device
    poscount: fiterwheel actual position count
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

### **def FWxC\_COMMAND\_LIB.FWxCGetSensorMode ( hdl, sensormode)**

```
get the fiterwheel current sensor mode
Args:
    hdl: the handle of opened FWxC device
    sensormode: fiterwheel actual sensor mode:0, Sensors turn off;1, Sensors remain active
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

### **def FWxC\_COMMAND\_LIB.FWxCGetSpeedMode ( hdl, spemode)**

```
get the fiterwheel current speed mode
Args:
    hdl: the handle of opened FWxC device
    spemode: 0,slow speed;1,high speed
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

### **def FWxC\_COMMAND\_LIB.FWxCGetTriggerMode ( hdl, triggermode)**

```
get the fiterwheel current position count
Args:
    hdl: the handle of opened FWxC device
    triggermode: fiterwheel actual trigger mode:0, input mode;1, output mode
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

### **def FWxC\_COMMAND\_LIB.FWxCIsOpen ( serialNo)**

```
Check opened status of FWxC device
Args:
    serialNo: serial number of FWxC device
Returns:
    0: FWxC device is not opened; 1: FWxC device is opened.
```

### **def FWxC\_COMMAND\_LIB.FWxCListDevices ()**

```
List all connected FWxC devices
Returns:
    The FWxC device list, each deice item is [serialNumber, FWxCType]
```

**def FWxC\_COMMAND\_LIB.FWxCOpen ( serialNo, nBaud, timeout)**

```
Open FWxC device
Args:
    serialNo: serial number of FWxC device
    nBaud: bit per second of port
    timeout: set timeout value in (s)
Returns:
    non-negative number: hdl number returned Successful; negative number: failed.
```

**def FWxC\_COMMAND\_LIB.FWxCSave ( hdl)**

```
save all the settings as default on power up
Args:
    hdl: the handle of opened FWxC device
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

**def FWxC\_COMMAND\_LIB.FWxCSetPosition ( hdl, pos)**

```
set fiterwheel's position
Args:
    hdl: the handle of opened FWxC device
    pos: fiterwheel position
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

**def FWxC\_COMMAND\_LIB.FWxCSetPositionCount ( hdl, count)**

```
set fiterwheel's position count
Args:
    hdl: the handle of opened FWxC device
    count: fiterwheel PositionCount
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

**def FWxC\_COMMAND\_LIB.FWxCSetSensorMode ( hdl, senmode)**

```
set fiterwheel's sensor mode
Args:
    hdl: the handle of opened FWxC device
    senmode: fiterwheel sensor mode
        sensors=0 Sensors turn off when wheel is idle to eliminate stray light
        sensors=1 Sensors remain active
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

**def FWxC\_COMMAND\_LIB.FWxCSetSpeedMode ( hdl, spmode)**

```
set fiterwheel's trigger mode
Args:
    hdl: the handle of opened FWxC device
    spmode: fiterwheel speed mode
        speed=0 Sets the move profile to slow speed
        speed=1 Sets the move profile to high speed
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```



**def FWxC\_COMMAND\_LIB.FWxCSetTriggerMode ( hdl, trimode)**


```
set fiterwheel's trigger mode
Args:
    hdl: the handle of opened FWxC device
    trimode: fiterwheel's trigger mode
        trig=0 Sets the external trigger to the input mode, Respond to an active low pulse by
        advancing position by 1
        trig=1 Sets the external trigger to the output mode, Generate an active high pulse when
        selected position arrived at
Returns:
    0: Success; 0xEA: CMD_NOT_DEFINED; 0xEB: time out; 0xED: invalid string buffer.
```

## Chapter 4 LabVIEW Software Development Kit

The user can start software development with LabVIEW 2013 or later versions based on LabVIEW instrument driver mechanism. The supported files are in *LabVIEW SDK* under the **Sample** directory.

### 4.1. How to install


Unzip the zip file and **copy** to instr.lib folder under LabVIEW installation folder.



Destination folder: under %LabVIEW install path%\instr.lib


Typically, C:\Program Files (x86)\National Instruments\LabVIEW 2013\instr.lib

Note: LabVIEW 2013 or later LabVIEW versions are compatible.



### 4.2. How to find VI


VI Could be found under: Functions\Instrument I/O\Instrument Drivers\



### 4.3. How to use

#### 1. From VI

**Note:** Before you open the SDK LabVIEW project, make sure the device has been connected to the computer.



#### 2. From VI tree

Some classic data flow in VI tree.

Use the Example Finder to find examples demonstrating the usage of this instrument driver.  
To launch Example Finder, select "Find Examples..." from the LabVIEW Help menu.

**EDIT:** Create example data file (.bin3) for Example Finder

#### Initialize



#### Open



#### Action/Status



#### Configuration



#### Data




#### Close



### 3. From example

An examples show the classic usage. Example path: instr.lib\Thorlabs\_FWxC\Examples



Easy programming and detailed comment will help.



**THORLABS**  
[www.thorlabs.com](http://www.thorlabs.com)

---